





Equipped with a head pillow and foldable side supports with safety lock, the ARJO Basic Shower Trolley is height adjustable, helping to reduce postural load.



The Shower Panel, available in three different sink alternatives, provides an efficiently controlled delivery system for assisted showering. The unit features an accurate thermostatic mixer with shut-off. As well as the resident shower, panels are fitted with an equipment cleaning shower to ensure good system hygiene routines.



ARJO, Inc. 50 N. Gary Avenue Roselle, Illinois 60172 Phone: 800-323-1245 Fax: 888-594-2756

www.arjo.com

The soft mattress is comfortable to lie on and the wide and sturdy construction of the trolley enables easy access when showering the resident. The trolley is easily adjusted from either side to a comfortable working height.



The ARJO Maxi Slide or Maxi Tube, when used in conjunction with the Basic Shower Trolley, simplifies the transfer from bed.

## **Product Information**

Height adjustable between 22.8" (570 mm) and 35.6" (890 mm)

All four castors have brakes

Equipped with head pillow and foldable side supports with safety lock

Length	76" (1900 mm)
Width	30.4" (760 mm)
Max load	330 lbs (150 kg)
Accesson	

Drain hose

## **Resident Gallery**

## Suitable for Emma



Emma, who is almost completely bedridden and totally dependant

Please contact ABJO for further information on the Resident Gallery<sup>16</sup>.

Use only ARUO Designed Parts, i.e. parts specifically designed for the purpose, on equipment and products supplied by ARUO. Our Conditions of Sale specify that ARUO is not liable for injuries attributable to the use of inadequate parts. Continuous development is part of our policy, and we therefore reserve the right to change specifications without notice.

Basic is a trademark belonging to ARJO International AG.



ISO 9001 • BS 5750 Part 1 • EN 29001 ABJO Hospital Equipment AB: Sweden

MEMBER OF THE GETINGE GROUP